

## SPECIFICATIONS

<b>Product name</b>	<b>SNEEZER BEAM</b>
Model	SN-211
Technology of treatment	Dual wavelength phototherapy
Illuminate position	Nostrils
Wavelengths used	657-940 nanometer
Emitting mode	Continuous mode
Battery type	9V, heavy duty dry battery (1 piece)
Battery life	Approximately 100 times (approximately 3 minutes each time)
Product dimension	78 x 65 x 31 mm
Product weight	80 gr.
Operating temperature/humidity	Temperature: 15~40 °C (59 °~104 °F); Humidity: 30% to 80%
Storage temperature/humidity	Temperature: -10~60 °C (14 °~140 °F) Humidity: 30% to 80%
Accessories	Main unit, 1 piece 9V battery, instruction manual
Warranty	24 months from the date of purchase
Certification	ISO-13485, EMC, ISO-10993, EN60601-1, CE 0434.