

SPECIFICATIONS

Product name	NEPLUS – SUPER MESH NEBULIZER
Model	NE-SM1
Power	Lithium ion battery 3,6~4,2V DC 1,4W
Dimensions	52 x 53 x 109 mm
Weight	Main unit 88 g.
Particle size	Average 5 micron
Nebulization rate	Approx. 0,5 ml/min (23°C)
Volume of medication bottle	6 ml
Operating time with full charge	4 hours
Charging life cycle of battery	Over 300 cycles
Operating environment	Temp. 10°C~40.0°C, Humidity 30%~85%
Storage environment	Temp. -20°C~60°C, Humidity 10%~95%
Warranty	2 years from purchase
Certification	CE0120, ISO, TUV